 **axoguard**
nerve connector

Distributed by:

 **axogen**

Axogen Corporation
13631 Progress Blvd., Suite 400
Alachua, Florida 32615 USA
1.888.Axogen1 (1.888.296.4361)
1.386.462.6800
customer care@axogeninc.com
www.axogeninc.eu

Manufactured by:
Cook Biotech Incorporated
1425 Innovation Place
West Lafayette, Indiana 47906 U.S.A.

©2021 Axogen Corporation

CE
2797

EC REP

EC Representative:
Cook Biotech Europe ApS
Sandet 6, DK-4632
Bjæverskov, DENMARK
Tel: +45 56868686

CH REP

Swiss Representative:
Decomplex AG
Freiburgstrasse 3
3010 Bern
Switzerland

October 2021

FP0131-01D

Symbols Used on Labeling

Caution
Atención
Vorsicht
Attention
Attenzione
Forsigtig
Försiktighet
Forsiktig
Let op
Προσοχή
Aviso
Dikkat
Figyelem
Pozor
Uwaga

Use-by date
Fecha de caducidad
Verfallsdatum
Utiliser avant
Data di scadenza
Anvendes inden
Används före
Utløpsdato
Uiterste gebruksdatum
Ημερομηνία λήξης
Data de validade
Son kullanna tarihi
Lejárati idő
Datum použitelnosti
Termin ważności

Do not re-use
No reutilizar
Nicht zur Wiederverwendung
Ne pas réutiliser
Non riutilizzare
Må ikke genanvendes
Får inte återanvändas
Skal ikke gjenbrukes
Niet opnieuw gebruiken
Μην επαναχρησιμοποιείτε
Não reutilizar
Tekrar kullanmayın
Tilos ismételtlen felhasználni
Nepoužívejte opakovaně
Nie używać powtórnie

Batch code
Código de lote
Chargenbezeichnung
Code de lot
Codice del lotto
Batchkode
Satsnummer
Kode for parti
Batchcode
Κωδικός παρτίδας
Código de lote
Parti Kodu
Gyártási tétel kódja
Kód šarže
Numer serii

Manufacturer
Fabricante
Hersteller
Fabricant
Fabbicante
Producent
Tillverkare
Produsent
Fabrikant
Κατασκευαστής
Fabricante
Üretici
Γυάρτό
Výrobce
Producent

STERILE EO

Sterilized using ethylene oxide
Esterilizado con óxido de etileno
Sterilisation mit Ethylenoxid
Stérilisé à l'oxyde d'éthylène
Sterilizzato con ossido di etilene
Steriliseret med ethylenoxid
Steriliserad med etylenoxid
Sterilisert med etylenoksid
Gesteriliseerd met ethyleenoxide
Αποστειρώθηκε με χρήση οξειδίου του αιθυλενίου
Esterilizado com óxido de etileno
Etilen oksitle sterilize edilmiştir
Etilén-oxiddal sterilizálva
Sterilizováno ethylenoxidem
Sterylizowany tlenkiem etylenu

Keep dry
Mantener seco
Vor Nässe schützen
Conserver au sec
Tenere al riparo dall'umidità
Opbevares tørt
Förvaras torrt
Holdes tørr
Droog houden
Διατηρείτε στεγνό
Manter seco
Kuru Tutun
Szárason tartandó
Chraňte před vlhkem
Chronić przed wilgocią

Magnetic resonance safe
Seguro en entornos de resonancia magnética (MR Safe)
Magnetresonanzsicher (MR Safe)
Sans danger avec la résonance magnétique (MR Safe)
Sicuro per la risonanza magnetica (MR Safe)
MR-sikker (MR Safe)
MR-säker (MR Safe)
Magnetresonanssikker (MR Safe)
Veilig voor MRI (MR Safe)
Ασφαλές σε περιβάλλον μαγνητικού συντονισμού (MR Safe)
Seguro para ressonância magnética (MR Safe)
Manyetik Rezonans Güvenli (MR Safe)
MR-biztonságos (MR Safe)
Bezpečný při vyšetření MRI (MR Safe)
Bezpieczny w środowisku rezonansu magnetycznego (MR Safe)

REF

Catalogue number
Número de catálogo
Bestellnummer
Numéro de catalogue
Numero di catalogo
Katalognummer
Katalognummer
Katalognummer
Catalogusnummer
Αριθμός καταλόγου
Número de catálogo
Katalog numarası
Katalógusszám
Katalogové číslo
Numer katalogowy

EC REP

European authorized representative
Representantes autorizados en Europa
Bevollmächtigter in der EU
Mandataire agréé dans l'Union européenne
Rappresentante autorizzato per l'Europa
Autoriseret repræsentant i EU
Auktoriserad representant i Europa
Autoriseret representant i Det europeiske fællesskab
Gemachtigde in Europa
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
Representante autorizado da Comunidade Europeia
Avrupa onaylı temsilci
Hivatalos európai képviselő
Zplnomocněný zástupce v Evropě
Autoryzowany przedstawiciel na terenie Europy

Date of manufacture
Fecha de fabricación
Herstellungsdatum
Date de fabrication
Data di fabbricazione
Fremstillingsdato
Tillverkningsdatum
Produksjonsdato
Fabricagedatum
Ημερομηνία κατασκευής
Data de fabricaçã
Üretim tarihi
Gyártás dátuma
Datum výroby
Data produkcji

CH **REP**

Swiss Authorized Representative
Representante autorizado en Suiza
Schweizer Bevollmächtigter
Mandataire sis en Suisse
Mandatario svizzero
Autoriseret repræsentant i Schweiz
Auktoriserad representant i Schweiz
Autorisert representant i Sveits
Zwitserse gemachtigde
Εξουσιοδοτημένος αντιπρόσωπος στην Ελλάδα
Representante autorizado na Suíça
İsviçre Yetkili Temsilcisi
Meghatalmazott képviselő Svájcban
Zplnomocněný zástupce ve Švýcarsku
Upoważniony przedstawiciel w Szwajcarii

Description

The Axoguard Nerve Connector is a surgical implant designed to align and connect peripheral nerves. Axoguard Nerve Connector is designed to be an interface between the nerve and the surrounding tissue, thus providing a protective environment for natural peripheral nerve healing. Axoguard Nerve Connector is comprised of an extracellular matrix (ECM) and is fully remodeled during the healing process. Animal model data suggests that the Axoguard Nerve Connector is quickly vascularized in situ and gradually remodeled into a connective tissue structure similar to nerve epineurium over several months¹. When hydrated, Axoguard Nerve Connector is easy to handle, soft, pliable, nonfriable, and porous. Axoguard Nerve Connector is flexible to accommodate movement of the joint and associated tendons, and has sufficient mechanical strength to hold sutures. Axoguard Nerve Connector is provided sterile, for single use only, and in a variety of sizes to meet the surgeon's needs.

Indications for Use

The Axoguard Nerve Connector is indicated for the repair of peripheral nerve discontinuities with gaps up to 5 mm. The Axoguard Nerve Connector is supplied sterile and is intended for single use.

Rx ONLY

This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

This product is intended for use by trained medical professionals.

Contraindications

The Axoguard Nerve Connector is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

NOTE: This device is not intended for use in vascular applications.

Precautions

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize device.**
- Discard all open and unused portions of the device.
- Device is sterile provided the package is dry, unopened and undamaged. Do not use device if the package seal is damaged or open.
- Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Do not suture device prior to rehydration.

Potential Complications

Possible complications can occur with any nerve repair surgical procedure including pain, infection, decreased or increased nerve sensitivity, and complications associated with use of anesthesia.

If any of the following conditions occur and cannot be resolved, careful removal of the device should be considered:

- Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction
- Infection

Storage

The Axoguard Nerve Connector should be stored in a clean, dry location at room temperature.

Sterilization

This device has been sterilized with ethylene oxide.

Suggested Instructions for Use

NOTE: Always handle Axoguard Nerve Connector using aseptic technique. Minimize contact with latex gloves.

1. Follow standard operating procedures for exposure and mobilization of the nerve (see Figure 1). Determine the nerve diameter in millimeters (mm) using a suitable measuring instrument. Select an Axoguard Nerve Connector of sufficient diameter to account for normal edema following traumatic nerve injury and to allow easy insertion of the nerve stumps into the device.

2. Open the outer carton and remove the sterile pouch. Using standard aseptic technique, open the pouch and pass the inner tray to the sterile field for further handling.
3. Open the tray and, if necessary, trim the Axoguard Nerve Connector to a length that is appropriate for the nerve gap. Allow for a maximum 5 mm nerve gap and at least 2 mm on the proximal and distal stumps for insertion into the connector. The Axoguard Nerve Connector may be trimmed after hydration if desired.
4. Fill the pre-molded rehydration reservoir with room temperature sterile saline or sterile Lactated Ringer's solution. Hydrate the Axoguard Nerve Connector for 10 seconds or until the desired handling characteristics are achieved, but not more than 20 minutes.
5. Hemostasis of both nerve stumps must be achieved prior to beginning the entubulation procedure. In the event a tourniquet is used, release the tourniquet and achieve hemostasis before beginning the entubulation procedure.
6. Use non-absorbable suture to secure the Axoguard Nerve Connector in place. Pass the suture through the wall of the tube from outside to inside, at least 1 mm from the edge (see Figure 2). Then pass the suture transversely through the epineurium of one nerve stump at a distance of at least 1 mm from the cut nerve face. Reverse the suture and pass it through the wall of the Axoguard Nerve Connector by pulling the suture so that the nerve stump is drawn into the tube. Tie the suture so that it is secure but avoid generating tension in the suture itself. If desired, add an additional suture to secure the device.

7. Use a syringe to gently flush the lumen of the Axoguard Nerve Connector with sterile saline or Lactated Ringer's solution (Figure 3). Repeat the suturing procedures described above with the opposite nerve stump (Figure 4). Fill the interior of the Axoguard Nerve Connector with saline or Lactated Ringer's solution (Figure 5). See Figure 6 for completed repair.

NOTE: Clinical literature has shown that better outcomes are found with using conduits with a gap of 5 mm or less.^{2,3}


Figure 3

Figure 4


Figure 5

Figure 6

8. Discard any unused portions of the Axoguard Nerve Connector according to institutional guidelines for biological waste. Do not resterilize.

How Supplied

Axoguard Nerve Connector is provided in a plastic tray and outer sterile pouch. The pouch is heat-sealed to provide a sterile barrier and has a peelable seal. Contents of the package are guaranteed sterile unless the package is opened or damaged. The Axoguard Nerve Connector and packaging do not contain natural rubber latex. *Do not use if the Axoguard Nerve Connector peel pouch appears to be open or damaged.*

References

1. Kokkalis ZT, Pu C, Small GA, Weiser RW, Venouziou AI, Sotereanos DG. Assessment of processed porcine extracellular matrix as a protective barrier in a rabbit nerve wrap model. *J Reconstr Microsurg.* 2011;27(1):19-28.
2. Isaacs J, Safa B, Evans PJ, Greenberg J. Technical assessment of connector assisted nerve repair. *J Hand Surg.* 2016;41(7):760-766.
3. Ducic I, Safa B, DeVinney E. Refinements of nerve repair with connector-assisted coaptation. *Microsurgery.* 2017;37(3):256-263.

Inquiries

For additional information, to place an order, or to report adverse events, contact:

Axogen Customer Care: 888-Axogen1 (888-296-4361)

E-mail: customer care@axogeninc.com

Returned Goods Policy

Authorization from Axogen Customer Care must be obtained prior to returning product. Sterile product must be returned in unopened, undamaged cartons, packed to prevent damage.

‡ Axogen Corp. owns registrations for, or other trademark rights in the “a” mark; AXOGEN; AXOGUARD; AXOGUARD NERVE CONNECTOR; and AXOGEN, REVOLUTIONIZING THE SCIENCE OF NERVE REPAIR in various countries throughout the world.

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Referencias

- Kokkalis ZT, Pu C, Small GA, Weiser RW, Venouziou AI, Sotereanos DG. Assessment of processed porcine extracellular matrix as a protective barrier in a rabbit nerve wrap model. *J Reconstr Microsurg.* 2011;27(1):19-28.
- Isaacs J, Safa B, Evans PJ, Greenberg J. Technical assessment of connector assisted nerve repair. *J Hand Surg.* 2016;41(7):760-766.
- Ducic I, Safa B, DeVinney E. Refinements of nerve repair with connector-assisted coaptation. *Microsurgery.* 2017;37(3):256-263.

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Referanser

- Kokkalis ZT, Pu C, Small GA, Weiser RW, Venouziou AI, Sotereanos DG. Assessment of processed porcine extracellular matrix as a protective barrier in a rabbit nerve wrap model. *J Reconstr Microsurg.* 2011;27(1):19-28.
- Isaacs J, Safa B, Evans PJ, Greenberg J. Technical assessment of connector assisted nerve repair. *J Hand Surg.* 2016;41(7):760-766.

3. Ducic I, Safa B, DeVinney E. Refinements of nerve repair with connector-assisted coaptation. *Microsurgery*. 2017;37(3):256-263.

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1. Kokkalis ZT, Pu C, Small GA, Weiser RW, Venouziou AI, Sotereanos DG. Assessment of processed porcine extracellular matrix as a protective barrier in a rabbit nerve wrap model. *J Reconstr Microsurg.* 2011;27(1):19-28.
2. Isaacs J, Safa B, Evans Pd, Greenberg J. Technical assessment of connector assisted nerve repair. *J Hand Surg.* 2016;41(7):760-766.
3. Ducic I, Safa B, DeVinney E. Refinements of nerve repair with connector-assisted coaptation. *Microsurgery.* 2017;37(3):256-263.

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axogen, revolutionizing the science of nerve repair †