# axoguard nerve protector

Distributed by:

 axogen

Axogen Corporation  
13631 Progress Blvd., Suite 400  
Alachua, Florida 32615 USA  
1.888.Axogen1 (1.888.296.4361)  
1.386.462.6800  
customer care@axogeninc.com  
www.axogeninc.eu

Manufactured by:  
Cook Biotech Incorporated  
1425 Innovation Place  
West Lafayette, Indiana 47906 U.S.A.

©2021 Axogen Corporation

**CE**  
2797

**EC REP**

EC Representative:  
Cook Biotech Europe ApS  
Sandet 6, DK-4632  
Bjæverskov, DENMARK  
Tel: +45 56868686

**CH REP**

Swiss Representative:  
Decomplex AG  
Freiburgstrasse 3  
3010 Bern  
Switzerland

## Symbols Used on Labeling

Caution  
Atención  
Vorsicht  
Attention  
Attenzione  
Forsigtig  
Försiktighet  
Forsiktig  
Let op  
Προσοχή  
Aviso  
Dikkat  
Figyelem  
Pozor  
Uwaga

Use-by date  
Fecha de caducidad  
Verfallsdatum  
Utiliser avant  
Data di scadenza  
Anvendes inden  
Används före  
Utløpsdato  
Uiterste bruksdatum  
Ημερομηνία λήξης  
Data de validade  
Son kulanma tarihi  
Lejáratí idő  
Datum použitelnosti  
Termin ważności

Do not re-use  
No reutilizar  
Nicht zur Wiederverwendung  
Ne pas réutiliser  
Non riutilizzare  
Må ikke genanvendes  
Får inte återanvändas  
Skal ikke gjenbrukes  
Niet opnieuw gebruiken  
Μην επαναχρησιμοποιείτε  
Não reutilizar  
Tekrar kullanmayın  
Tilos ismételtlen felhasználni  
Nepoužívejte opakovaně  
Nie używać powtórnie

Batch code  
Código de lote  
Chargenbezeichnung  
Code de lot  
Codice del lotto  
Batchkode  
Satsnummer  
Kode for parti  
Batchcode  
Κωδικός παρτίδας  
Código de lote  
Parti Kodu  
Gyártási tétel kódja  
Kód šarže  
Numer serii

Manufacturer  
Fabricante  
Hersteller  
Fabricant  
Fabbricante  
Producent  
Tillverkare  
Produsent  
Fabrikant  
Κατασκευαστής  
Fabricante  
Üretici  
Gyártó  
Výrobce  
Producent

**STERILE EO**

Sterilized using ethylene oxide  
Esterilizado con óxido de etileno  
Sterilisation mit Ethylenoxid  
Stérilisé à l'oxyde d'éthylène  
Sterilizzato con ossido di etilene  
Steriliseret med ethylenoxid  
Steriliserad med etylenoxid  
Sterilisert med etylenoksid  
Gesteriliseerd met ethyleenoxide  
Αποστειρώθηκε με χρήση οξειδίου του αιθυλενίου  
Esterilizado com óxido de etileno  
Etilen oksitle sterilize edilmiştir  
Etilén-oxidál sterilizálva  
Sterilizováno ethylenoxidem  
Sterylizowany tlenkiem etylenu

Keep dry  
Mantener seco  
Vor Nässe schützen  
Conserver au sec  
Tenere al riparo dall'umidità  
Opbevares tørt  
Förvaras torrt  
Holdes tørr  
Droog houden  
Διατηρείτε στεγνό  
Manter seco  
Kuru Tutun  
Szárason tartandó  
Chraňte před vlhkem  
Chronić przed wilgocią

Magnetic resonance safe  
Seguro en entornos de resonancia magnética (MR Safe)  
Magnetresonanzsicher (MR Safe)  
Sans danger avec la résonance magnétique (MR Safe)  
Sicuro per la risonanza magnetica (MR Safe)  
MR-sikker (MR Safe)  
MR-säker (MR Safe)  
Magnetresonanssikker (MR Safe)  
Veilig voor MRI (MR Safe)  
Ασφαλές σε περιβάλλον μαγνητικού συντονισμού (MR Safe)  
Seguro para ressonância magnética (MR Safe)  
Manyetik Rezonans Güvenli (MR Safe)  
MR-biztonságos (MR Safe)  
Bezpečný při vyšetření MRI (MR Safe)  
Bezpieczny w środowisku rezonansu magnetycznego (MR Safe)

**REF**

Catalogue number  
Número de catálogo  
Bestellnummer  
Numéro de catalogue  
Numero di catalogo  
Katalognummer  
Katalognummer  
Katalognummer  
Catalogusnummer  
Αριθμός καταλόγου  
Número de catálogo  
Katalog numarası  
Katalógusszám  
Katalogové číslo  
Numer katalogowy

**EC REP**

European authorized representative  
Representantes autorizados en Europa  
Bevollmächtigter in der EU  
Mandataire agréé dans l'Union européenne  
Rappresentante autorizzato per l'Europa  
Autoriseret repræsentant i EU  
Auktoriserad representant i Europa  
Autorisert representant i Det europeiske fellesskap  
Gemachtigde in Europa  
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα  
Representante autorizado da Comunidade Europeia  
Avrupa onaylı temsilci  
Hivatalos európai képviselő  
Zplnomocněný zástupce v Evropě  
Autoryzowany przedstawiciel na terenie Europy

Date of manufacture  
Fecha de fabricación  
Herstellungsdatum  
Date de fabrication  
Data di fabbricazione  
Fremstillingsdato  
Tillverkningsdatum  
Produksjonsdato  
Fabricagedatum  
Ημερομηνία κατασκευής  
Data de fabricaçoã  
Ûretim tarihi  
Gyártás dátuma  
Datum výroby  
Data produkcji

CH REP

Swiss Authorized Representative  
Representante autorizado en Suiza  
Schweizer Bevollmächtigter  
Mandataire sis en Suisse  
Mandatario svizzero  
Autoriseret repræsentant i Schweiz  
Auktoriserad representant i Schweiz  
Autorisert representant i Sveits  
Zwitserse gemachtigde  
Εξουσιοδοτημένος αντιπρόσωπος στην Ελλάδα  
Representante autorizado na Suíça  
İsviçre Yetkili Temsilcisi  
Meghatalmazott képviselő Svájcban  
Zplnomocněný zástupce ve Švýcarsku  
Upoważniony przedstawiciel w Szwajcarii

**Description**

The Axoguard Nerve Protector is a surgical implant that provides non-constricting protection for peripheral nerves. Axoguard Nerve Protector is designed to be an interface between the nerve and the surrounding tissue, thus providing a protective environment for natural peripheral nerve healing. Axoguard Nerve Protector is comprised of an extracellular matrix (ECM) and is fully remodeled during the healing process. Animal model data suggests that the Axoguard Nerve Protector is quickly vascularized in situ and gradually remodeled into a connective tissue structure similar to nerve epineurium over several months<sup>1</sup>. When hydrated, Axoguard Nerve Protector is easy to handle, soft, pliable, nonfriable, and porous. Axoguard Nerve Protector is flexible to accommodate movement of the joint and associated tendons, and has sufficient mechanical strength to hold sutures. Axoguard Nerve Protector is provided sterile, for single use only, and in a variety of sizes to meet the surgeon's needs.

**Indications for Use**

The Axoguard Nerve Protector is indicated for the repair of peripheral nerve injuries where there is no gap. The Axoguard Nerve Protector is supplied sterile and is intended for single use.

**Rx ONLY**

This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

This product is intended for use by trained medical professionals.

**Contraindications**

The Axoguard Nerve Protector is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

*NOTE: This device is not intended for use in vascular applications.*

**Precautions**

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize device.**
- Discard all open and unused portions of the device.
- Device is sterile provided the package is dry, unopened and undamaged. Do not use device if the package seal is damaged or open.
- Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Do not suture device prior to rehydration.

**Potential Complications**

Possible complications can occur with any nerve repair surgical procedure including pain, infection, decreased or increased nerve sensitivity, and complications associated with use of anesthesia.

If any of the following conditions occur and cannot be resolved, careful removal of the device should be considered:

- Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction
- Infection

## Storage

The Axoguard Nerve Protector should be stored in a clean, dry location at room temperature.

## Sterilization


This device has been sterilized with ethylene oxide.

## Suggested Instructions for Use


*NOTE: Always handle Axoguard Nerve Protector using aseptic technique. Minimize contact with latex gloves.*

1. Follow standard operating procedures for exposure and mobilization of the nerve. Determine the nerve diameter in millimeters (mm) using a suitable measuring instrument. If necessary, repair the nerve using standard operating procedures. Select an Axoguard Nerve Protector of sufficient diameter to account for normal edema following traumatic nerve injury. The Axoguard Nerve Protector diameter should be at least 1-2 mm larger than the measured nerve diameter and long enough to cover the affected area.
2. Open the outer carton and remove the sterile pouch. Using standard aseptic technique, open the pouch and pass the inner tray to the sterile field for further handling.
3. Open the tray and, if necessary, trim the Axoguard Nerve Protector to the appropriate dimensions for covering the damaged portion of the nerve. The Axoguard Nerve Protector may be trimmed after hydration, if desired.
4. Fill the pre-molded rehydration reservoir with room temperature sterile saline or sterile Lactated Ringer's solution. Hydrate the Axoguard Nerve Protector for 10 seconds or until the desired handling characteristics are achieved, but not more than 20 minutes.
5. Position the Axoguard Nerve Protector around the nerve (Figure 1). If transected, re-approximate the nerve ends (Figure 2). If desired, gently flush the device with sterile saline or Lactated Ringer's solution to improve conformability with the nerve. Secure the device as necessary. This may include placing running sutures along the longitudinal slit to enclose the nerve, and/or stay sutures through the nerve epineurium away from the coaptation site. See Figure 3 for completed repair.

### Repair of peripheral nerve injuries where there is no gap

### Re-approximate nerve ends

6. Discard any unused portions of the Axoguard Nerve Protector according to institutional guidelines for biological waste. Do not resterilize.

## How Supplied

Axoguard Nerve Protector is provided in a plastic tray and outer sterile pouch. The pouch is heat-sealed to provide a sterile barrier and has a peelable seal. Contents of the package are guaranteed sterile unless the package is opened or damaged. The Axoguard Nerve Protector and packaging do not contain natural rubber latex. Do not use if the Axoguard Nerve Protector peel pouch appears to be open or damaged.

## References

1. Kokkalis ZT, Pu C, Small GA, Weiser RW, Venouziou AI, Sotereanos DG. Assessment of processed porcine extracellular matrix as a protective barrier in a rabbit nerve wrap model. *J Reconstr Microsurg.* 2011;27(1):19-28.

## Inquiries

For additional information, to place an order, or to report adverse events, contact:

**Axogen Customer Care: 888-Axogen1 (888-296-4361)**

**E-mail: [customercare@axogeninc.com](mailto:customercare@axogeninc.com)**

## Returned Goods Policy

Authorization from Axogen Customer Care must be obtained prior to returning device. Sterile product must be returned in unopened, undamaged cartons, packed to prevent damage.

‡ Axogen Corp. owns registrations for, or other trademark rights in the “a” mark; AXOGEN; AXOGUARD; AXOGUARD NERVE PROTECTOR; and AXOGEN, REVOLUTIONIZING THE SCIENCE OF NERVE REPAIR in various countries throughout the world.

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