ENGLISH

## Instructions for Use

**DESCRIPTION**  
Avance® Nerve Graft is processed nerve allograft intended for bridging nerve discontinuities in peripheral nerve regeneration. It is desiccated and cleaned extracellular matrix from donated human peripheral nerve. The cleaning process preserves the inherent microstructural structure of the native nerve extracellular matrix, while clearing the cells, cellular debris and certain proteins sufficient to cytotoxicity (CSPC).

Avance® Nerve Graft is implanted and connects the proximal and distal ends of a transected nerve. The processed nerve allograft serves as a scaffold allowing regenerating axons to grow into the patient's distal nerve tissue towards the target muscle or skin.

Avance® Nerve Graft offers structural characteristics and handling similar to nerve autograft; pliability of soft tissue, an intact epineurium to suture the nerve graft in place, and intact endoneurial tubes for the axons to grow through.

Avance® Nerve Graft is supplied sterile in a variety of lengths and diameters to allow the surgeon to tailor the size of the grafted nerve and provide a tensionless repair of the defect. It is for single patient use only. Approximate graft lengths and diameters are listed on the package label.

**REGULATORY CLASSIFICATION**  
Avance® Nerve Graft is processed and distributed in accordance with US FDA requirements for Human Cellular and Tissue-based Products (HCT/P) under 21 CFR Part 1271 regulations. US State regulators and the standards of the American Association of Tissue Banks (AATB). Additionally, international regulations are followed as appropriate.

Avance® Nerve Graft is to be dispensed only by or on the order of a licensed physician.

**INDICATIONS FOR USE**  
Avance® Nerve Graft is processed nerve allograft (human) intended for the surgical repair of peripheral nerve discontinuities to support regeneration across the defect.

**CONTRAINDICATIONS**  
Avance® Nerve Graft is contraindicated for use in any patient in whom soft tissue implants are contraindicated. This includes any pathway that would limit the blood supply and compromise healing or evidence of a current infection.

**WARNINGS**  
Caution against screening laboratory testing, tissue processing, and gamma irradiation have been used to minimize the risks of transmission of relevant communicable diseases to the patient. As with any processed human donor tissue, Avance® Nerve Graft cannot be guaranteed to be free of all pathogens.

Do not reuse or re-sterilize Avance® Nerve Graft and do not refreeze the graft once it has been thawed.

**DONOR RECOVERY AND SCREENING**  
Each donor consent for donation of peripheral nerve tissue for the recovery of the peripheral nerve tissue is performed in an aseptic manner by FDA registered and state licensed (where required) US tissue banks. Donor eligibility is carefully evaluated as required by the US FDA and US State regulators. Additionally, donor eligibility is determined in accordance with AATB standards and appropriate international regulations. Tissue donors are evaluated for high risk behaviors and communicable diseases as determined by the Centers for Medicine and Medical Services (CMS). Testing is conducted using test kits approved by the US FDA. As disease screening methods for certain diseases may not be detected. The following complications are:

- Transmission of diseases of unknown origin;
- Transmission of known infectious agents including, but not limited to viruses, bacteria, and fungi;

**INCQUIRIES**  
For additional information, to place an order, or to report errors, accidents or adverse events contact:

Avance® Customer Care  
Phone: 888.AxoGen1 (888.296.4361) or 386.462.8800  
Email: CustomerCare@AxoGen.com

Customers outside of US: Contact the AxoGen authorized distributor servicing your facility, email AxoGen® Customer Care or contact AxoGen directly in the US at 386.462.8800.

Additional testing may be performed, as required by local authorities in international markets. If required, each donor is tested and shown to be **negative or nonreactive** for the following:

- Human T-Cell Lymphotropic Virus (HTLV) Type I Antibody
- Human T-Gall Lymphotropic Virus (HTLV) Type II Antibody

All testing is performed by a laboratory registered with FDA to perform donor testing and testing of human tissue. Testing is performed by the US Clinical Laboratory Improvement Amendments (CLIA) of 1988 and 42 CFR Part 493, or that has met equivalent requirements as determined by the Centers for Medicine and Medical Services (CMS). The equivalent requirements are followed in accordance with ISO 15171 guidelines.

The Medical Director of AxoGen (US state licensed) has determined that the tissue is suitable for transplantation. Records of all testing and medical releases are maintained by AxoGen.

**PROCESSING**  
Avance® Nerve Graft is processed in controlled environments using Good Tissue Practice procedures to prevent contamination and ensure maximum containment of the tissue. Processing involves the use of proprietary physiological buffers, enzymes and surfactants and the processed tissue may contain traces of these processing agents. The cleaning process preserves the inherent microstructural structure of the native nerve extracellular matrix while clearing the cells, cellular debris and certain proteins such as chondroitin sulfate proteoglycans (CSPC). After completion of processing, Avance® Nerve Graft is sized, packaged and sterilized using gamma irradiation in accordance with ISO 15171 guidelines.

Avance® Nerve Graft has been tested in accordance with ISO 10993 standards. The test results demonstrated that the processed nerve allograft is biocompatible.

The processed nerve allograft tissue will naturally vary in color from white, off-white, pink, pale pink, and yellow to pale yellow. Occasional dark spots or localized discoloration is a normal occurrence.

**HOW SUPPLIED**  
Avance® Nerve Graft is placed into a plastic tray and then inserted into chevron pouches. Each chevron pouch is heat-sealed to provide a moisture barrier and prevent a chevron seal. The outer chevron pouch is foil to provide a moisture barrier. Approximate graft lengths and diameters are listed on the package label. Avance® Nerve Graft is irradiated and supplied frozen. Contents are sterile until the package is opened or damaged.

**TRANSPORT AND STORAGE**  
Avance® Nerve Graft is shipped frozen from dry ice via validated shipping containers. Immediately upon receipt, verify presence of dry ice in shipper and place product in freezer at -40°C (-40°F). Keep the product frozen until the time of use. Avance® Nerve Graft is a registered trademark of AxoGen. AxoGen is a registered trademark of E.DuPont De Nemours & Co.

Avance® Nerve Graft is a registered trademark of AxoGen. AxoGen is a registered trademark of E.DuPont De Nemours & Co.

Avance® Nerve Graft has been tested in accordance with ISO 10993 standards. The test results demonstrated that the processed nerve allograft is biocompatible.

The processed nerve allograft tissue will naturally vary in color from white, off-white, pink, pale pink, and yellow to pale yellow. Occasional dark spots or localized discoloration is a normal occurrence.

**INSTRUCTIONS FOR USE**  
1. Follow standard operating procedures for exposure and mobilization of the injured peripheral nerve.

2. Determine the injured nerve diameter and length in millimeters (mm) using a suitable instrument. Measure the nerve diameter at the site of the injury.

3. Select Avance® Nerve Graft(s) of comparable diameter to match the native nerve and of sufficient length to ensure a tension free repair.

4. Remove the foil outer chevron pouch containing Avance® Nerve Graft. Instructions for Use, patient record labels and Avance® Nerve Graft (TUR) from the package.

5. Compare the diameter and length of the native nerve with the diameter and length on the package. If the numbers do not match, **DO NOT USE** the product and notify AxoGen Customer Care immediately. See product label for expiration date.

6. Using standard aseptic technique, peel open the outer foil chevron pouch and pass the inner Tyevek® chevron pouch to the sterile field for further handling.

7. Open the Tyevek® chevron pouch and remove the plastic tray.

8. Open the plastic tray and fill the pre-moistened thawing reservoir with room temperature sterile saline or sterile Lactated Ringers solution (LRS). Do not heat the graft or add heated saline to the reservoir.

9. Allow Avance® Nerve Graft to thaw completely before using which will take approximately 5 - 10 minutes. Once thawed, Avance® Nerve Graft is soft and pliable throughout. Avance® Nerve Graft is comprised of a sachet and outer sleeve of queleu maniere que ce soit. **NE PAS UTILISER** la greffe nerveuse Avance® Nerve Graft et informer immédiatement le service clientel d'AxoGen.

10. Complete and send the Tissue Utilization Report (TUR) back to AxoGen.

TISSUE UTILIZATION REPORT (TUR)  
Each Avance® Nerve Graft package contains a Tissue Utilization Report (TUR). In accordance with US FDA, US Joint Commission and international requirements, a TUR should be completed for each Avance® Nerve Graft used in the procedure and returned to AxoGen or other representative as described on the TUR.

Record the distinct HCT/P identification code in hospital or facility records and in the patient's file. Complete all information on the card, affix ONE (1) peel-off label of each Avance® Nerve Graft used, seal and return to AxoGen or other representative as described on the TUR.

**POTENTIAL COMPLICATIONS**  
The following adverse events are anticipated and may potentially occur during treatment:

- mild incisional redness;
- tenderness of surgical area;
- mild edema of surgical area;
- controllable pain at surgical area;
- decreased sensation; and,
- numbness

These adverse events are considered expected and are not required to be recorded unless they increase in severity.

Inherent risks of any surgical procedure include, infection, blood loss, and anesthesia associated complications. Complications specific to nerve reconstruction include pain, decreased or increased sensitivity, and impaired motor or sensory function. As with all peripheral nerve surgery, there is a risk of failure of the nerve to regenerate.

Hypersensitivity, allergic reactions, or other adverse immune responses have not been observed in preclinical studies or reported clinically in the use of Avance® Nerve Graft. Because Avance® Nerve Graft is composed of proteins such as collagen and laminin, there is potential may exist for such reactions.

Avance® Nerve Graft is processed human nerve tissue. As with all donated human tissue products, the risk for transmission of communicable disease does exist. Robust donor screening and selection criteria, completed as required by the FDA and in accordance with AATB, state, and federal guidelines, processing controls, and terminal sterilization with gamma irradiation greatly reduce but cannot totally eliminate this risk. Per AATB standards, communicable disease testing is performed by a laboratory registered with the FDA to perform donor testing and certified to perform such testing on human specimens under the US Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR part 493, or that has met equivalent requirements as determined by the Centers for Medicine and Medical Services (CMS). Testing is conducted using test kits approved by the US FDA. As disease screening methods for certain diseases may not be detected. The following complications may occur:

- Transmission of diseases of unknown origin;
- Transmission of known infectious agents including, but not limited to viruses, bacteria, and fungi;

**Single European Code (SEC) For Tissues**  
If applicable SEC Code is found on or attached to the outer carton.

**DISPOSAL**  
Disposal of Avance® Nerve Graft in accordance with local, state and federal or country regulations for disposal of human tissue.

**REFERENCES**  
1. ISO 11317:2006 Sterilization of health care products — Radiation guidelines  
2. ISO 10993:2002 Biological evaluation of medical devices

**REF**  
Product Code

**LOT**  
Lot Number

**TEMPERATURE LIMIT**  
Do Not Reuse/Restilize

**STERILE R**  
Contents are sterile unless outer package is damaged, Radiation Sterilization

**EXP**  
Expiration Date

**REF**  
Product Code

**LOT**  
Lot Number

**TEMPERATURE LIMIT**  
Do Not Reuse/Restilize

**STERILE R**  
Contents are sterile unless outer package is damaged, Radiation Sterilization

**EXP**  
Expiration Date

**REF**  
Product Code

**LOT**  
Lot Number

**TEMPERATURE LIMIT**  
Do Not Reuse/Restilize

**STERILE R**  
Contents are sterile unless outer package is damaged, Radiation Sterilization

**EXP**  
Expiration Date

**REF**  
Product Code

**LOT**  
Lot Number

**TEMPERATURE LIMIT**  
Do Not Reuse/Restilize

**STERILE R**  
Contents are sterile unless outer package is damaged, Radiation Sterilization

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